



PubMed Nucleotide Protein Genome Structure PopSet Taxonomy OMIM Books

Search PubMed for Go Clear

Limits Preview/Index History Clipboard Details

Display Abstract Sort Save Text Add to Clipboard Order

1: Biotechniques 2000 Jan;28(1):107-109

Related Articles, Books

High-throughput method for isolating plasmid DNA with reduced lipopolysaccharide content.

Neudecker F, Grimm S.

Max-Planck-Institute for Biochemistry, Martinsried, Germany.

Isolating plasmid DNA from bacteria is a fundamental step in molecular biology. It is often accomplished by an alkaline lysis of bacteria and the subsequent adsorption of nucleic acids to silica oxide in the presence of chaotropic substances. Here we show that the addition of such chaotropic reagents is not required for the efficient DNA isolation with silica oxide. This surprising finding allowed us to purify plasmid DNA with significantly less lipopolysaccharides (LPS), which is otherwise a common bacterial contaminant of silica oxide-isolated DNA and inhibits subsequent applications. In addition, we have implemented a precipitation step that altogether leads to a reduction of the LPS content by a factor of 900 relative to published methods. Our novel protocol facilitates an inexpensive high-throughput analysis of pure plasmids in a 96-well format without the addition of hazardous reagents.

Publication Types:

- Technical report

PMID: 10649779 [PubMed - indexed for MEDLINE]

Display Abstract Sort Save Text Add to Clipboard Order

[Write to the Help Desk](#)
[NCBI](#) | [NLM](#) | [NIH](#)
Department of Health & Human Services
[Freedom of Information Act](#) | [Disclaimer](#)